






# Round Clamp (WC-R1) or (WC-R7/8)

Scale 1 : 1



## Set Includes:

- (2) Clamps
- (2) 1/4-28 x 2 1/4" Socket Head Cap Screws
- (2) 1/4-28 Lock-nuts



**Mount'n Mover™**  
by BlueSky Designs

[mountnmover.com](http://mountnmover.com)

Product:  
Weight:

WC-R1" or WC-R7/8" Round Clamp  
5.6 oz. (0.16 kg)

1-888-724-7002 (toll free) or 1-612-724-7002

# Wheelchair Attachment Overview

\* All plates are considered 'multi-use plates' and can be configured, combined, and used in a variety of ways.



Start with a power or manual wheelchair

## Anchor to Chair

\* What you use will depend on what kind of chair frame you have.

- Round tubing
- Slide track
- Hole-spacing
- Tie-Down



T-Nut (WC-TN)



Adapter Plate 2 (WC-AP2)



Adapter Plate 4 (WC-AP4)




Round Clamp (WC-R1 or R7/8)

Span the bridge clamps with a plate




Bridge Clamps



Adapter Plate 2 (WC-AP2)




(WC-AP)



(WC-AP4)

## OPTIONAL OFFSET



Adapter Plate 2 (WC-AP2)



Adapter Plate 4 (WC-AP4)



Angle Adjust Plate (WC-APAP)



L-Angle Extension Plate (WC-LAE)

## Wheelchair Bracket



(WB2)

Or



(WB)

\* Offset example





Plate to span across the clamps